



# HYGRO-GUARD / CLIMA- / DATA-LOG

Professional, accurate long run data recording - fixed and mobile

## HUMIDITY- / TEMPERATURE- / PRESSURE-MEASUREMENT

### QUICK FACTS

## HYGRO-GUARD/CLIMA- / DATA-LOG

High capacity data logging systems for data mapping, monitoring and validation

Data logger series with 3.2 million data point memory. Delivered ready for use including batteries, factory calibration certificate, data visualisation and programming software and data cable. Available in various models and with different humidity measurement technologies.

### ADVANTAGES & BENEFIT

## BIG DATA STORAGE AND DISPLAY, INTUITIVE VISUALISATION SOFTWARE, DATA EXPORT, NETWORK-COMPATIBLE

Autonomous battery powered solution, which records data for a very long period. Very flexible, ideal for network integration or for stand-alone use. Robust instrument housing for harsh and mobile applications. Available also with high precision electrolytic humidity sensor (DataLog).

### FEATURES

	HygroGuard	ClimaLog	DataLog
BATTERY POWERED	X	X	X
ALTERNATIVELY POWERED WITH MAINS SUPPLY	X	X	X
NETWORK-COMPATIBILITY	X	X	X
INTEGRATED HUMIDITY & TEMPERATURE SENSORS	X	X	
INTEGRATED ABSOLUTE PRESSURE SENSOR		X	
CONNECTABLE EXTERNAL SENSORS			X
CAPACITIVE HUMIDITY MEASUREMENT	X	X	
ELECTROLYTIC HUMIDITY MEASUREMENT (EXT.)			X
3.2 MILLION DATA POINTS STORAGE CAPACITY	X	X	X
INCL. VISUALISATION DATA EXPORT SOFTWARE	X	X	X
LARGE LC-DISPLAY	X	X	X

### AREA OF APPLICATIONS

ROOM MAPPING IN PHARMACEUTICAL INDUSTRY  
 LOGISTICS AND WAREHOUSES  
 CALIBRATION LABS  
 CLIMATIC CHAMBERS  
 BUILDING DRYING  
 ARCHIVES, MUSEUMS, RIPENING CELLARS  
 AND MORE...

**novasina**  
 The Art of Precision Measurement